


An SP Industries Company

Lab Liquid Mixing Equipment

• LG-MS-PRO Magnetic Stirrers with Heating



RTD Temperature Probe with ± 0.2 °C accuracy

RS232 Port

Low Profile Design

Large LCD Display

LG-MS-PRO

LG-MS-PRO Technical Data

Number of Stirring Positions	1
Stirring Bar Length Max (mm)	80
Stirring Quantity Max (L)	20
Speed Range (rpm)	100-1500
Heat Output (W)	500
Adjustable Temperature Max (°C)	340
Control Accuracy with Sensor (°C)	± 0.2
RS232 Serial Port	Yes
External Temperature Sensor	Yes
Dimensions [DxWxH] (mm)	280x160x85



Silicone Protective Cover



Optional Quarter System for Reaction Vials

Other Models

Product Number	Description	Price
LG-MS-H-PRO	Magnetic Hot Plate Stirrer, Stainless Steel Top, w/o Temperature Probe	\$455.00
LG-MS-H-PRO-C	Magnetic Hot Plate Stirrer, Ceramic Top, w/o Temperature Probe	\$455.00



LG-MS-H-PRO-C

Product Number	Included Items	Price
LG-MS-PRO	<ul style="list-style-type: none"> ✓ Magnetic Stirrer Hot Plate, Stainless Steel Top ✓ PT1000 RTD Temperature Probe with Stand ✓ Silicone Protective Cover 	<div style="border: 2px solid black; border-radius: 50%; padding: 5px; display: inline-block;"> ONLY \$601.52 </div>



Toll Free: **1-800-220-5171**
Intl.: **856-691-3200**



Toll Free: **1-800-220-1081**
Intl.: **856-691-6206**



www.wilmad-labglass.com
cs@wilmad-labglass.com

Temperature Probe

Product Number	Description	Price
LG-PT-1000	RTD Temperature Probe	\$125.00
LG-SUPPORT	LG-PT-1000 Stand Set	\$58.50




Recommended for use with LG-MS-H-PRO and LG-MS-H-PRO-C

Quarter System for Vial

Product Number	Description	Depth (mm)	# of Vials	Color	Price
LG-1890001	Ø135 mm Carrier Plate			Blue	\$67.45
LG-1890003	Quarter, 20 mL Reaction Vessel	24	4	Purple	\$43.26
LG-1890004	Quarter, 30 mL Reaction Vessel	30	4	Blue	\$47.38
LG-1890005	Quarter, 40 mL Reaction Vessel	43	Black	\$51.09	
LG-1890002	Quarter, 4 mL Reaction Vessel	19	11	Red	\$53.05
LG-1890048	Quarter, 8 mL Reaction Vessel	26	6	Green	\$51.24
LG-1890049	Quarter, 16 mL Reaction Vessel	32	4	Golden	\$47.38



LG-1890003




LG-1890004




LG-1890005



LG-1890002



LG-1890048



LG-1890049

Recommended for use with stainless steel top hot plates, which include LG-MS-PRO and LG-MS-H-PRO

Reaction Block for Flask

Product Number	Description	# of flasks	Price
LG-1890080	Reaction Block for 50 mL Round Bottom Flask	1	\$119.00
LG-1890081	Reaction Block for 100 mL Round Bottom Flask	1	\$124.00
LG-1890082	Reaction Block for 250 mL Round Bottom Flask	1	\$136.00
LG-1890083	Reaction Block for 500 mL Round Bottom Flask	1	\$143.00



LG-1890080



LG-1890081



LG-1890082



LG-1890083

Recommended for use with stainless steel top hot plates, which include LG-MS-PRO and LG-MS-H-PRO