

Lab Liquid Mixing Equipment

• LG-MX Vortex Mixer



LG-MX Technical Data	
Power (W)	60
Shaking Movement	orbital
Motor Input (W)	58
Motor Output (W)	10
Speed Range (rpm)	0-2500
Switch	Continuous/Touch
Dimension [D×W×H] (mm)	127×130×160

- Easily Set Speed from 0 to 2500 rpm
- Select from Continuous On/Off/Touch On
- Water-proof Switch Cover
- Standard Anticorrosive Silicone Base
- Exchangeable Vacuum Feet Included for Improved Stability

Product Number	Included Items	Price
LG-MX	<ul style="list-style-type: none"> ✓ Vortex Mixer with Standard Silicone Base ✓ Exchangeable Vacuum Feet (4×) 	



Toll Free: **1-800-220-5171**
Intl.: **856-691-3200**



Toll Free: **1-800-220-1081**
Intl.: **856-691-6206**



www.wilmad-labglass.com
cs@wilmad-labglass.com