



NEW
Wilmad Sets the Standard

Low Pressure/Vacuum Shigemi Tube for Bruker, Agilent (Varian) and JEOL Spectrometers

SP Service VirTis Hull Wilmad | LabGlass Genevac FTS Hotpack



See back for various solvent-matched models, accessories and pricing
Shigemi is a trademark of Shigemi, Inc.



ISO 9001:2008 Certified
An SP Industries Company



Phone



Fax



Web

Toll Free: **1-800-220-5171**
Intl.: **856-691-3200**

Toll Free: **1-800-220-1081**
Intl.: **856-691-6206**

www.wilmad-labglass.com
cs@wilmad-labglass.com

5 mm O.D. Standard and Gas-Tight Shigemi Tube

Product No.	Description	O.D. (mm)	Bottom Length (mm)	Matched Solvent	Compatibility	Price Each
CMS-005B	Standard Shigemi tube-full set ¹	4.965	8	Chloroform-d	Bruker	\$98.44
CMS-005J		4.965	12		JEOL	\$98.44
CMS-005V		4.965	15		Agilent(Varian)	\$98.44
CMS-005B-LPV	Gas-tight Shigemi tube-full set ²	4.965	8		Bruker	\$241.77
CMS-005J-LPV		4.965	12		JEOL	\$241.77
CMS-005V-LPV		4.965	15		Agilent(Varian)	\$241.77
MMS-005B	Standard Shigemi tube-full set ¹	4.965	8	Methanol-d ₄	Bruker	\$110.75
MMS-005J		4.965	12		JEOL	\$110.75
MMS-005V		4.965	15		Agilent(Varian)	\$110.75
MMS-005B-LPV	Gas-tight Shigemi tube-full set ²	4.965	8		Bruker	\$259.96
MMS-005J-LPV		4.965	12		JEOL	\$259.96
MMS-005V-LPV		4.965	15		Agilent(Varian)	\$259.96
DMS-005B	Standard Shigemi tube-full set ¹	4.965	8	DMSO-d ₆	Bruker	\$110.75
DMS-005J		4.965	12		JEOL	\$110.75
DMS-005V		4.965	15		Agilent(Varian)	\$110.75
DMS-005B-LPV	Gas-tight Shigemi tube-full set ²	4.965	8		Bruker	\$259.96
DMS-005J-LPV		4.965	12		JEOL	\$259.96
DMS-005V-LPV		4.965	15		Agilent(Varian)	\$259.96
BMS-005B	Standard Shigemi tube-full set ¹	4.965	8	Deuterium Oxide	Bruker	\$98.44
BMS-005J		4.965	12		JEOL	\$98.44
BMS-005V		4.965	15		Agilent(Varian)	\$98.44
BMS-005B-LPV	Gas-tight Shigemi tube-full set ²	4.965	8		Bruker	\$241.77
BMS-005J-LPV		4.965	12		JEOL	\$241.77
BMS-005V-LPV		4.965	15		Agilent(Varian)	\$241.77
529-C	Positioning Rod			Universal	All	\$19.47
529-B	PTFE Holder				All	\$28.63

- Standard Shigemi tube full set includes one outer tube and one inner insert
- Gas-tight Shigemi tube full set includes one outer tube fused with complete LPV valve, one inner insert with detachable PTFE Holder (529-B) and one Positioning Rod (529-C)

10 mm O.D. Standard Shigemi Tube¹

Product No.	Description	O.D. (mm)	Bottom Length (mm)	Matched Solvent	Compatibility	Price Each
CMS-010B	Standard Shigemi tube-full set ²	9.98	8	Chloroform-d	Bruker	\$332.35
CMS-010V	Standard Shigemi tube-full set ²	9.98	15		Agilent(Varian)	\$363.40
MMS-010B	Standard Shigemi tube-full set ²	9.98	8	Methanol-d ₄	Bruker	\$345.00
MMS-010V	Standard Shigemi tube-full set ²	9.98	15		Agilent(Varian)	\$374.90
DMS-010B	Standard Shigemi tube-full set ²	9.98	8	DMSO-d ₆	Bruker	\$345.00
DMS-010V	Standard Shigemi tube-full set ²	9.98	15		Agilent(Varian)	\$374.90
BMS-010B	Standard Shigemi tube-full set ²	9.98	8	Deuterium Oxide	Bruker	\$332.35
BMS-010V	Standard Shigemi tube-full set ²	9.98	15		Agilent(Varian)	\$363.40

- Please contact us for the pricing on 10 mm O.D. Gas-tight Shigemi tube
- Standard Shigemi tube full set includes one outer tube and one inner insert