

Committed to Product Excellence and Continuous Support

an SP Industries Company

Wilmad-LabGlass Equipment



CE ISO 9001/13485
CERTIFIED

Full Line of Liquid Transfer Devices

- Flowing Smoothly at Precision Volume

1172 NW Boulevard, Vineland, NJ 08360, USA
Tel: 800-220-5171 Fax: 800-220-1081
www.wilmad-labglass.com

Med Plus Pipettor



- The permissible errors of Med Plus Pipettor are on average 2 times better than those recommended by ISO.
- Each Med Plus is individually tested and comes with a certificate.

Product No.	Volume Range	Volume (μ l)	Med Plus Maximum Permissible Errors		ISO 8655 Maximum Permissible Errors		Price Each
			Inaccuracy	Imprecision	Inaccuracy	Imprecision	
LG-MPP2.5	0.1-2.5 μ l	0.25	\pm 12.00%	\leq 6.00%	\pm 40%	\leq 20%	\$92
		1.25	\pm 3.00%	\leq 3.00%	\pm 8%	\leq 4%	
		2.5	\pm 2.50%	\leq 2.00%	\pm 4%	\leq 2%	
LG-MPP10	0.5-10 μ l	1	\pm 2.50%	\leq 1.50%	\pm 12%	\leq 8%	\$92
		5	\pm 1.50%	\leq 1.50%	\pm 2.4%	\leq 1.6%	
		10	\pm 1.00%	\leq 0.80%	\pm 1.2%	\leq 0.8%	
LG-MPP20	2-20 μ l	2	\pm 3.00%	\leq 2.00%	\pm 10%	\leq 5%	\$92
		10	\pm 1.2%	\leq 1.00%	\pm 2%	\leq 1%	
		20	\pm 0.90%	\leq 0.40%	\pm 1%	\leq 0.5%	
LG-MPP100	10-100 μ l	10	\pm 3.00%	\leq 1.50%	\pm 8%	\leq 3%	\$92
		50	\pm 1.00%	\leq 0.40%	\pm 1.6%	\leq 0.6%	
		100	\pm 0.80%	\leq 0.15%	\pm 0.8%	\leq 0.3%	
LG-MPP200	20-200 μ l	20	\pm 3.00%	\leq 1.00%	\pm 8%	\leq 3%	\$92
		100	\pm 0.80%	\leq 0.30%	\pm 1.6%	\leq 0.6%	
		200	\pm 0.60%	\leq 0.15%	\pm 0.8%	\leq 0.3%	
LG-MPP1000	100-1000 μ l	100	\pm 0.80%	\leq 0.30%	\pm 1.6%	\leq 0.6%	\$92
		500	\pm 0.70%	\leq 0.25%	\pm 1.6%	\leq 0.6%	
		1000	\pm 0.60%	\leq 0.20%	\pm 0.8%	\leq 0.3%	

MedSpenser Plus Dispenser

- Same precision as Eppendorf's Varispenser/ Varispenser Plus.
- Each MedSpenser Plus is individually tested and comes with a certificate.

Med Plus Maximum Permissible Errors

Product No.	Volume Range	Graduation (ml)	Inaccuracy	Imprecision	Price Each
LG-MPD5	0.5-5 ml	0.1	±0.5%	0.1%	\$154
LG-MPD10	5-50 ml	1.0	±0.5%	0.1%	\$231

Each dispenser comes with 25, 28, 32, 38 and 40 mm thread adapters



MedFiller Plus Electronic Pipette Filler




Product Number	LG-MFP
Price	\$195
Aspiration Speed	>5 ml/sec
Dispense Speed	>5 ml/sec to gravity dispense
Battery	Lithium-Ion
Battery Life	>8 hours when fully charged
Charging Time	~2 hours
Pipette Types	Glass or Plastic (0.1-100 ml), Pasteur pipettes
Filter	Standard and sterile - 0.20 µm and 0.45 µm

Manual Pipette Controller

- Reduce hand stress compared to conventional squeeze bulb.

Product Number	Color	Price
LG-MCP-B	Blue	\$31
LG-MCP-R	Red	\$31
LG-MCP-G	Green	\$31
LG-MCP-Y	Yellow	\$31
LG-MCP-P	Pink	\$31

See back cover for color options



Manual Pipette Controller



an SP Industries Company

Committed to Product Excellence and Continuous Support

SP Service

VirTis

Hull

Wilmad | LabGlass

Genevac

FTS

Hotpack

SP Industries Corporate Headquarters:

935 Mearns Road Warminster, PA 18974 USA Tel: 800.523.2327 215.672.7800 Fax: 215.672.7807 www.spindustries.com