

Section 8


PRECISION BORE SQUARE TUBING

Wilmad square precision bore tubing is characterized by the same accuracy and uniformity that marks all Wilmad products. The sizes listed are available from tooling on hand in our plant.

The dimensions shown in the tables are for Pyrex Brand, Duran or KG-33 Glass. The same tooling can be used to produce tubes from different types of glass but there will be a slight change in dimensions.

SQUARE PB TUBING

Dimensions	Length inches	Length mm
8mm x 8mm	11	279.40
10mm x 10mm	10	254.00
13mm x 13mm	12	304.80
16mm x 16mm	12	304.80
20mm x 20mm	12	304.80
1 inch x 1 inch	14	355.60
1 1/8 in. x 1 1/8 in.	12	304.80
1 1/4 in. x 1 1/4 in.	12	304.80
1 1/2 in. x 1 1/2 in.	12	304.80



We can also produce tubing in shapes other than square to precision bore tolerances. Such shapes include rectangular, triangular and other geometric configurations. For your special requirements, call, email or Fax our Engineering Department with details.