



1172 NW Boulevard Vineland NJ 08360

800.220.5171

856.691.3200

Fax: 800.220.1081

856.691.6206

www.wilma-labglass.com

Technical Information: Round Bottom Flask Diameters

Capacity (mL)	Diameter (mm)	Capacity (mL)	Diameter (mm)
5	26	250	85
10	31	300	87
15	36	500	105
20	39	1,000	131
25	41	2,000	166
35	46	3,000	185
50	49	5,000	223
65	55	12,000	295
100	62	22,000	345
125	69	50,000	460
200	76	72,000	520